

Simultaneous Twist Decapper for 96 multi tubes!

Simultaneous De-capping & Re-capping for 96(48) Screw caps

© CREATE DREAM

- ★ High Speed :16 seconds for 96(48) decapping: 20 seconds for recapping
- ★ Applicable for SBS format 96(48) well screw caps and tubes
- ★ Compatible with Pippeting device for integration (TD-96(48)R)
- ★ Adjustable for Manual access (TD-96(48),TD-96R(48)) & Robotic transport (TD-96(48)R)

Two types are available!

Stand Alone Type:TD-96 (48) with internal controller

Integration Type:TD-96 (48) R with external controller

Adjustable for various types of transporting & pipetting systems



Stand Alone Type
for manual operation
TD-96 (48)



Integration Type
for pipetting devices
TD-96R (48)



Manufacturer :

Distributor :



IS Technology Japan, Inc.

5-9-1 Tokodai, Tsukuba, Ibaraki, 300-2635 Japan

TEL:+81-29-848-3308 FAX:+81-29-848-3309

Web site :<http://www.istj.co.jp>


Simultaneous Twist Decapper for 96 multi tubes!

© CREATE DREAM

Specifications

Product code	TD-96 (48)	TD-96 (48) R
Operation Type	Manual	Automation
Tube rack Insertion	Push button · Start	External controller
Decap & Recap	Push button · Start	External controller
Tube rack ejection	Push button · Start	External controller
Controller	Internal	External
Dimensions	W=350,D=290,H=450mm、26kg	W=300,D=220,H=450mm、35kg (Excluding Transport part & controller)
Power source	AC 100V 4A (can change the voltage)	AC 100V 4A (Excluding Transport part & controller)

Optional Item (TD-96)



Special spacer: Applicable for short and middle size tubes

Manufacturer :

Distributor :



IS Technology Japan, Inc.

5-9-1 Tokodai, Tsukuba, Ibaraki, 300-2635 Japan

TEL:+81-29-848-3308 FAX:+81-29-848-3309

Web site :<http://www.istj.co.jp>