


PCR Plate Positioner 「P3」



For 384well PCR Plate

Compatible with SBS format and every manufacturer's PCR plates

Size

■ Height	14.0mm	■ Outer	anodized aluminum Coating
■ Width	127.7mm	■ Material	Aluminum
■ Depth	85.5mm	■ Weight	150g

Application

- Automated liquid handling systems using 384 plates has been difficult because the rigidity of 384 well plates were low. However, our product, Positioner, enables researchers to use 384 plates for the automated systems by increasing the rigidity and keeping them horizontal.
- With Positioner, easy sealing for PCR Plate can be achieved with a sealer because it can fix plates firmly.
- Positioner is also applicable for Shaker for micro well plates.
- Positioner is adjustable and customized for your PCR Plates.

Manufacturer :

IS Technology Japan, Inc.

5-9-1 Tokodai, Tsukuba, Ibaraki, 300-2635 Japan




TEL:+81-29-848-3308 FAX:+81-29-848-3309

Web site :<http://www.istj.co.jp>

Distributor :



PCR Plate Positioner 「P3」

Compatible Plates for 384well PCR Plate Positioner				
Positioner品番		P3-100(No rib)	P3-200(with small ribs)	P3-300(with ribs)
Positioner図番		J07-0668-M01002-1	J08-0739-M01001	J09-0870-M07001
COLOR		White	Blue	Black
Price(1ea)		¥23,000	¥42,000	¥50,000
Photo				
Plate Maker	PlateNo.			
日本ジエネティクス (Sorenson)	38971	×	△	△
日本ジエネティクス (Sorenson)	39691	○	×	○
日本ジエネティクス (ABgene)	AB-1384	×	○	N/A
日本ジエネティクス (ABgene)	TF-0384	○	○	N/A
日本ジエネティクス (ABgene)	AB-0937	N/A	○	N/A
サーモフィッシャー (ナルジエヌク / ABgene)	TF-0384	○	○	N/A
グライナー (Griner)	785201	×	○	N/A
グライナー (Griner)	785290	△	×	N/A
アプライドバイオシステムズ (ABI)	4309849	○	○	N/A
4titude		N/A	N/A	N/A
AxyGen		N/A	N/A	N/A

○:Perfect fitting , ×:No fitting (lose) , △:Good fitting(obstinate) , N/A:Not tested



Manufacturer :

IS Technology Japan, Inc.

5-9-1 Tokodai, Tsukuba, Ibaraki, 300-2635 Japan

TEL:+81-29-848-3308 FAX:+81-29-848-3309

Web site :<http://www.istj.co.jp>

Distributor :