





# PLATE SPINNER



## iSPIN-04-28

Front View



Top View



- Put SBS standard plate to 60mm in height to four.
- Can spin-down even unbalance of 100g.
- Perform number of revolutions and the acceleration by the control of the PC.
- By the adoption of the AC servomotor, the precision of the stop position is less than  $\pm 0.1\text{mm}$ .

Model No.	iSPIN-04-28S-RSIF	iSPIN-04-28SC-RSIF	iSPIN-00-RA01-RSIF
Type/Style	Standard	Cooling Unit	Robot Arm Option
Revolutions	2800 rpm	2800 rpm	-
IBT	100g	50g	-
Accuracy	$\pm 0.1\text{mm}$	$\pm 0.1\text{mm}$	$\pm 0.1\text{mm}$
Operation	Host Communication	Host Communication	Host Communication
Controller	Distinction	Distinction	Distinction
Dimensions	W=500, D=500, H=640mm, 145kg (Controller does not include)	W=500, D=500, H=640mm, 145kg (Controller does not include)	W=500, D=300, H=300mm, 16kg
Interlocks	Rack detection	Rack detection	Rack detection
Utility	AC 100V 15A	AC 100V 20A	AC 100V 5A